THUASNE

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KNEE IMMOBILISATION SPLINT - 3 PANELS

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

The device is composed of :

- 4 removable, radiolucent, rigid stays,
- 3 panels, adjustable to various leg circumferences,
- 7 self-fastening closure and support straps,
- Foam padding for patient comfort.

Bilateral product.

Available in 4 sizes.

Composition

Textile components: Polyester - Polyurethane - Cotton - Polypropylene - Polyamide.

Rigid components: Aluminium - Polyamide.

Properties/Mode of action

Knee joint stabilisation via the rigid stays.

Total immobilisation at 0°.

Adjustable closure and support via self-fastening straps.

Indications

Post-operative immobilisation.

Post-traumatic immobilisation.

Knee pain management.

Contraindications

Do not use the product if the diagnosis has not been confirmed. Do not use in the event of known allergy to any of the components. Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.

History of venous or lymphatic disorders.

Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

Do not use in the event of unstable fractures, fractures of the lower part of the tibia or fibula, or fractures of the upper part of the femur.

Precautions

Verify the product’s integrity before every use.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

It is recommended that a healthcare professional supervises the first application.

Strictly comply with your healthcare professional’s prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

It is recommended that the venous thromboembolic risk be assessed before any immobilisation is indicated. Follow the advice of a healthcare professional.

For hygiene and performance reasons, do not re-use the product for another patient.

It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.

In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.

In the event of any modification in the product’s performance, remove it and consult a healthcare professional. Prolonged immobilisation can sometimes cause muscle weakness. Do not wear the product when driving a vehicle.

Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

Undesirable side-effect

This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity. Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional supervises the first application.

It is recommended to wear the device directly on the skin, unless contraindicated.

- Open the splint and lay it flat open 1.

- Place the opened splint under the leg 2, the wider edge towards the top of the leg.

The patella must be centred in the cut-away area.

N.B. : The top of the splint must not impede flexion of the leg at the groin.

- Shaping the stays: if needed, it is possible to shape the stays in order to avoid potential pressure points or the splint slipping. Apply a slight pressure in the desired direction.

Repeat the process on the other stays if needed. If subsequent shaping of the stays is needed, consult a healthcare professional.

- Position the side panels in order to fully cover the leg (panels are adjustable to the thigh and calf circumference) 3.

- Engage the self-fastening straps in the buckles. Tighten the straps in the order described (A, B, C, D, E, F, G) 4. Once all the straps are fastened, readjust the tightening. If needed, reduce the strap length by cutting the edges with scissors.

Care/Maintenance

Product can be washed in accordance with the conditions shown on this instruction leaflet and on the label. Hand wash. Close the self-fastening tabs before washing. Remove stays before washing. Stays can be washed with a moist cloth. Do not use detergents,

fabric softeners or aggressive products (products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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© Thuasne - 6164801 (2020-12)